




# 1 billion frames per second and beyond





Delivering individual, specific  
ultra high-speed solutions that  
meet our clients needs exactly,  
whatever the application.

**That's our commitment.**



## Our company


Specialised Imaging is an internationally renowned company that focuses on the design and manufacture of ultra-high-speed imaging cameras for industrial, scientific and defence research applications.

The company was formed in 2003, its founder members having previously worked together in the high-speed imaging field and bringing over 80 years' combined experience to the venture.

Specialised Imaging has successfully launched many new and innovative ultra-high-speed imaging systems.

The company is at the forefront of world-wide innovation in the high-speed imaging field, having won the BEEA's Small Company of the Year award in 2009, the Queen's Award for Enterprise in 2011 and 2016.


This commitment to development has enabled the company to establish a reputation as an exciting and creative player in the high-speed camera market.





Innovative imaging solutions  
that incorporate the latest  
technological advances.

**That's our passion.**





## **Supporting you... and your camera**

At Specialised Imaging we relish new technological challenges, and we enjoy creating effective solutions. Producing a system that exactly meets your requirements, demands a company prepared and able to create specific optimised solutions.

Specialised Imaging has a strong track record in working with clients to design and develop new functions and facilities that fulfil their requirements.


This level of commitment and support continues throughout the life of your product – on-going advice, problem-solving and the design and reconfiguration of software are all part of our after-sales service.

## **World-wide partnerships**

To provide you with a total ultra-high-speed imaging solution, we have formed strong, strategic relationships with manufacturers in a variety of related fields. These strategic partnerships enable us to offer fully optimised imaging systems that include illumination, optical components, supports and triggering devices.

Close cooperation with our major suppliers also allows us to access and incorporate custom-designed components to specifically enhance the performance of our systems

Our reputation today extends worldwide with customers in North America, Asia-Pacific, Europe and the Middle East regions.



## Duplex multi channel framing camera

**Up to 1 Billion frames  
per second capture speed**

**36lp/mm system  
resolution**

**1360 x 1024 pixel,  
12-bit sensor  
resolution**

**Up to 16 discrete  
intensified optical  
channels**



The Specialised Imaging SIMD Framing Camera offers up to 32 images without creating shading, or parallax. Highly accurate timing and fully flexible intensified CCD sensors provide almost infinite control over interframe time, gain and exposure to capture even the most difficult ultra-fast phenomena.

Comprehensive triggering adjustment and a wide range of output signals are controlled using the custom software package which also includes measurement and image enhancement functions.

The SIMD has an optional port for the addition of a high-speed video camera to allow longer duration and simultaneous image capture. The Duplex camera configuration allows the number of images captured to be twice the number of channels.

### FEATURES

- ☐ Fully adjustable interframe time to 1ns
- ☐ Fully adjustable exposure down to 3ns
- ☐ Gain adjustment up to 7,000X
- ☐ Adjustable output triggers
- ☐ Nikon lens mount fitting
- ☐ Ethernet communications
- ☐ Duplex configuration camera

## MODELS

	<b>SIMD4</b>	<b>SIMD8</b>	<b>SIMD16</b>	<b>SIMD24</b>	<b>SIMD32</b>
Number of Channels	2	4	8	12	16
Number of images	4	8	16	24	32

## OPTICAL

Optics	Single input beam splitting optics Channels can be fitted with individual filters
Lenses	Nikon F-Mount
Internal electro-mechanical iris	f2.8 - f22
Shutter	Electro-mechanical
Distortion	Nominally zero
Channel Registration	Within one pixel after software correction
Intensity Variation	Better than 5% across the image
Auxiliary Optical Channel Interface	Nikon F-mount bayonet (Optional)

## INTENSIFIER / SENSOR

Image Sensor	ICX285AL (Intensified)
Active CCD Pixel	1360 (H) x 1024 (V)
Pixel Size	6.45 µm (H) x 6.45 µm (V)
Dynamic Range	12 bits
Intensifier	18mm High resolution MCP Input window Fused Silica Output window Fibre Optic Photocathode S25, others on request Phosphor screen P46
Gain	Variable up to 7,000
System resolution	>36 lp/mm

## MECHANICAL

Dimension cm (w/d/h)	22.5 x 63.0 x 53.0 (8CH, without lens)
Mount	3/8-16 UNC Female
Weight	24Kg (8CH without lens)

## TIMING PARAMETERS

System Clock	1GHz quartz crystal controlled
Exposure Mode (each image)	Single exposure or multiple exposures (Max. 8) per channel
Exposure Time	3ns - 10ms in 1ns steps independently variable
Separation Time (multiple exposure mode)	30ns - 20ms in 1ns steps independently variable
Interframe Time	0ns - 20ms in 1ns steps independently variable
Delay to 1st exposure	65ns to 10ms in 1ns steps, independently variable
Flash Outputs	5ns - 1ms in 1ns steps independently variable
Framing rates	up to 1 Billion fps

## INPUT / OUTPUT SIGNALS

Trigger 1	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Trigger 2	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Timing Monitor Pulses	Pulse width (min. 3ns) and position user programmable TTL into 50Ω
Flash Trigger Outputs	Pulse width (min. 5ns) and position user programmable TTL into 50Ω
Camera control	Data and command transfer via 100Mbps ethernet cable length 10m (standard), other lengths up to 100m
Software	Custom software compatible with Microsoft Windows Operating Systems for camera control, image data archiving in various file formats.
Power Requirements	100-240V AC 2A, 50-60Hz

## ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10 - 90% RH non condensing
Vibration shock	10 - 40 Hz Max. 10g in any direction
EMC	Meets all EC harmonized standards

**UK (Head Office / Factory)**  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com)

[info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## High resolution multi channel framing camera

**Up to 1 Billion frames per second capture speed**

**50lp/mm system resolution**

**1360 x 1024 pixel, 12-bit sensor resolution**

**Up to 16 discrete intensified optical channels**



The Specialised Imaging SIMX Framing Camera offers up to 16 high resolution images without creating shading, or parallax. Highly accurate timing and fully flexible intensified CCD sensors provide almost infinite control over interframe time, gain and exposure to capture even the most difficult ultra-fast phenomena.

Comprehensive triggering adjustment and a wide range of output signals are controlled using the custom software package which also includes measurement and image enhancement functions.

The SIMX has an optional port for the addition of a high-speed video, or streak camera to allow either simultaneous long duration or ultra high temporal resolution capture. A multi-spectral configuration SIMX camera can provide up to 16 different multi-spectral images with 5 colour and 1 monochrome images.

### FEATURES

- ☐ Fully adjustable interframe time to 1ns
- ☐ Fully adjustable exposure down to 3ns
- ☐ Gain adjustment up to 10,000X
- ☐ Adjustable output triggers
- ☐ Nikon lens mount fitting
- ☐ Ethernet communications
- ☐ Multi-Spectral configuration camera option



## MODELS

	<b>SIMX4</b>	<b>SIMX8</b>	<b>SIMX16</b>
Number of Channels	4	8	16

## OPTICAL

Optics	Single input beam splitting optics Channels can be fitted with individual filters
Lenses	Nikon F-Mount
Internal electro-mechanical iris	f2.8 - f22
Shutter	Electro-mechanical
Distortion	Nominally zero
Channel Registration	Within one pixel after software correction
Intensity Variation	Better than 5% across the image
Auxiliary Optical Channel Interface	Nikon F-mount bayonet (Optional)

## INTENSIFIER / SENSOR

Image Sensor	ICX285AL (Intensified)
Active CCD Pixel	1360 (H) x 1024 (V)
Pixel Size	6.45 $\mu$ m (H) x 6.45 $\mu$ m (V)
Dynamic Range	12 bits
Intensifier	18mm High resolution MCP Input window Fused Silica Output window Fibre Optic Photocathode S25, others on request Phosphor screen P43
Gain	Variable up to 10,000
System resolution	50 lp/mm

## MECHANICAL

Dimension cm (w/d/h)	22.5 x 63.0 x 53.0 (8CH, without lens)
Mount	3/8-16 UNC Female
Weight	24Kg (8CH without lens)

## TIMING PARAMETERS

System Clock	1GHz quartz crystal controlled
Exposure Mode (each image)	Single exposure or multiple exposures (Max. 8) per channel
Exposure Time	3ns - 10ms in 1ns steps independently variable
Interframe Time	0ns - 20ms in 1ns steps independently variable
Delay to 1st exposure	65ns to 10ms in 1ns steps, independently variable
Flash Outputs	5ns - 1ms in 1ns steps independently variable
Framing rates	up to 1 Billion fps

## INPUT / OUTPUT SIGNALS

Trigger 1	Electrical signal (BNC connector) Threshold variable from $\pm$ 25V Positive or Negative polarity, Make/Break 50 $\Omega$ or 1K $\Omega$ termination
Trigger 2	Electrical signal (BNC connector) Threshold variable from $\pm$ 25V Positive or Negative polarity, Make/Break 50 $\Omega$ or 1K $\Omega$ termination
Timing Monitor Pulses	Pulse width (min. 3ns) and position user programmable TTL into 50 $\Omega$
Flash Trigger Outputs	Pulse width (min. 5ns) and position user programmable TTL into 50 $\Omega$
Camera control	Data and command transfer via 100Mbps ethernet cable length 10m (standard), other lengths up to 100m available 100FX fibre optic ethernet link (up to 2Km) - optional
Software	Custom software compatible with Microsoft Windows Operating Systems for camera control, image data archiving in various file formats.
Power Requirements	100-240V AC 2A, 50-60Hz

## ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10 - 90% RH non condensing
Vibration shock	10 - 40 Hz Max. 10g in any direction
EMC	Meets all EC harmonized standards

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## The Ultimate High-Speed Video Camera

**Up to 5 Million images/  
second capture speed**

**924 x 768 pixel,  
10-bit sensor resolution**

**180 images**

**100ns minimum  
shutter**



The Specialised Imaging KIRANA is a true Ultra high-speed video camera that combines the flexibility of a video camera with the speed/resolutions approaching those only available with Framing cameras.

The unique custom design sensor offers 180 images at capture speeds up to 5 Million Images/second at full resolution.

In line with high-speed video cameras the KIRANA can be Frame synchronised with an external device such as another KIRANA or laser.

The KIRANA can also be recording prior to the event and triggered before, during or after the event.

### FEATURES

- ☐ Up to 5 Million images/second
- ☐ Adjustable exposure down to 100ns
- ☐ Pre & Post event triggering
- ☐ External synchronisation
- ☐ Nikon lens mount fitting
- ☐ Gigabit ethernet communications
- ☐ Compact and rugged design

**Kirana1M** Up to 1 Mfps

**Kirana5M** Up to 5 Mfps

### MODEL SPECIFICATION

	<b>Kirana1M</b>	<b>Kirana5M</b>
Frame Rate (Frames per second)	Up to 1 Mfps	Up to 5 Mfps
Exposure Time (minimum)	1µs 10ns step	100ns 10ns step
Trigger Mode	Start, End and Center	

### OPTICAL

Lenses	Nikon F-Mount
Shutter	Electro-mechanical
Distortion	Zero

### INTENSIFIER / SENSOR

Sensor	µCMOS (Non-intensified)
Number of Active Pixels	924 (W) x 768 (H)
Pixel Size	30µm
Digitisation	10bits Monochrome
Number of Frames	180

### MECHANICAL

Dimension mm (w/d/h)	<b>Head:</b> 22.8cm x 42cm x 19cm (without lens) <b>Power supply:</b> 19.5cm x 39.5cm x 19.5cm (inc. handle)
Weights	<b>Head:</b> 10.6Kg (23lbs) without lens. <b>Power Supply:</b> 4.8Kg (10.5lbs)
Head Mounting	3/8-16 UNC Female in head base.

### TIMING PARAMETERS

System Clock	200MHz quartz crystal controlled
Exposure time	100ns - 1/Frame rate
Framing rates	1000fps - 1Mfps or 5Mfps

### INPUT / OUTPUT SIGNALS

Trigger (2 off)	Electrical signal (BNC connector) Threshold variable from ± 25V Maximum Input level 50V Integral Velocity measurement system Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Video Out	XVGA
Aux Out	FSync or user programmable pulse width and position for strobe/laser illumination sources. TTL into 50Ω
Sync In	Input to allow the synchronisation of multiple cameras in Master-Slave configuration
Camera Control	Remote control via Standard 1Gbps Ethernet
Software	Custom software compatible with Microsoft Windows Operating Systems for control and data archiving in various file formats
Power Requirements	100-240V AC 2A, 50-60Hz
Saved Image Format	TIFF, JPEG, AVI or RAW

### ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10—90% RH non condensing
Vibration shock	10—40 Hz Max. 10g in any direction
EMC	Meets all EC harmonised standards


**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## Comprehensive High-Speed Linear Tracker System



**Award winning flight follower system**

**Multiple tracking modes**

**Remote control motorised adjustment**

**Multiple high-speed  
camera options**

The Specialised Imaging Tracker<sup>2</sup> is the next generation of projectile tracking platforms for high-speed video and measurement.

Full motorised remote control of three axis rotation and multiple inputs for real-time velocity adjustment contribute to the evolution of this award-winning system.

Built on a sturdy mount, the fully weatherproofed mirror and camera housings allow a large range of high-speed video cameras and long focal length lens options.

Custom software controls the Tracker system and provides calculators for Tracker placement, camera fields-of-view and velocities.

### FEATURES

- ☐ Full remote control operation
- ☐ Multiple operating modes allow capture of decelerating, accelerating, user defined velocity profiled projectiles
- ☐ Scan ratio range from 0.1 to 100
- ☐ Scanning accuracy of  $\pm 0.2^\circ$
- ☐ Gigabit ethernet communications
- ☐ Built in camera power, communications and trigger
- ☐ No calibration required

### OPERATING MODES

Fixed Velocity	Single trigger using known velocity
Velocity	The scan is corrected using the measured velocity from at least 2 of the 8 available detector inputs.
Position	The scan position is corrected from the detector inputs. Known velocity is assumed.
Drag / Acceleration	The scan is corrected using the measured velocity and drag / acceleration from at least 3 of the 8 available detector inputs.
Pre-defined profile	Programmable Velocity Vs Time curve. Triggered using single trigger. Used for non-linear projectile trajectories.
Advanced User Functions	Specialised Imaging is prepared to customise modes of operation to user requirements.
Skewed Geometry	Allows non perpendicular operation

### OPERATING PARAMETERS

Scan Ratio (SR)	0.1 to 100 (defined as the ratio of projectile velocity/ stand-off distance)
Scanning range (Max.)	-60° to +60 °
Scanning Distance	>=2x standoff distance (distance from the line of flight to Tracker2)
Scanning Accuracy	±0.2°
Positional Accuracy	±0.018°
Calibration	Not required
Projectile Velocity	SR x Standoff distance
Projectile Drag	0 to 100 m/s/m
Acceleration Angle	1° - 5° depending on scan rate (defined as the angle required to accelerate the mirror from rest to full scanning speed)

### ENVIRONMENTAL

Storage temperature	-10°C to +74°C
Operating temperature	-5°C to +50°C
Warmup Period	Not Required
Humidity	10 - 90% RH non-condensing
Operational vibration	10G, 10-40Hz Max, any direction
EMC	Meets all EC harmonized standards

### INPUT / OUTPUT SIGNALS

Detector In	BNC
Number of inputs	8
Trigger In	Rising or Falling Edge pulse Make/break
Camera Trigger	TTL positive pulse
Communication Interface	Data and command transfer via 1Gbps ethernet cable length 100m (standard). Other lengths available 1000FX fibre optic ethernet link (up to 2Km) - optional
Software	Custom software compatible with Microsoft Windows Operating Systems for control and data archiving in various file formats

### MECHANICAL

Dimensions mm (w/d/h)	1340 x 670 x 590 (without tripod)
Mount	Tripod Included

### MIRROR

Type	Optical flat elliptical Silicon Carbide Mirror
Size (HxW) mm	135 x 85 x 2

### CONTROL UNIT


System Clock	10MHz quartz crystal controlled
Trigger Jitter	<1us

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)





## Compact High-Speed Linear Tracker System



**Built on award winning flight  
follower system**

**Multiple tracking modes**

**Manual positional adjustment**

**Integrated high-speed  
camera options**

The Specialised Imaging Tracker<sup>Lite</sup> is a smaller, lighter system that retains the core functions of the Award winning Tracker<sup>2</sup> system.

Manual adjustment of three axis rotation and a single input for real-time velocity adjustment contribute to the simplification of the system.

Built on a standard tripod, the fully weatherproofed mirror and camera housing allow a dual line of high-speed video cameras and long focal length lens options.

Custom software controls the Tracker system and provides calculators for Tracker placement, camera fields-of-view and velocities.

### FEATURES

- ☐ Multiple operating modes allow capture of decelerating, accelerating, user defined velocity profiled projectiles
- ☐ Scan ratio range from 1 to 40
- ☐ Scanning accuracy of  $\pm 0.2^\circ$ .
- ☐ Gigabit ethernet communications
- ☐ Built in camera power, communications and trigger
- ☐ No calibration required

### OPERATING MODES

Fixed Velocity	Single trigger using known velocity
Velocity	The scan is corrected using the measured velocity from at least 2 of the 3 available detector inputs.
Position	The scan position is corrected from the detector inputs. Known velocity is assumed.
Drag	The scan is corrected using the measured velocity and drag from 3 detector inputs.
Pre-defined profile	Programmable Velocity Vs Time curve. Triggered using single trigger. Used for non-linear projectile trajectories.
Advanced User Functions	Specialised Imaging is prepared to customise modes of operation to user requirements.
Skewed Geometry	Allows non perpendicular operation

### OPERATING PARAMETERS

Scan Ratio (SR)	1 to 40 (defined as the ratio of projectile velocity/stand-off distance)
Scanning range (Max.)	-50° to +50°
Scanning Distance	>=2x standoff distance (distance from the line of flight to Tracker2)
Scanning Accuracy	±0.2°
Calibration	Not required
Projectile Velocity	SR x Standoff distance
Projectile Drag	0 to 100 m/s/m
Acceleration Angle	1° - 5° depending on scan rate (defined as the angle required to accelerate the mirror from rest to full scanning speed)

### ENVIRONMENTAL

Storage temperature	-10°C to +74°C
Operating temperature	-5°C to +45°C
Warmup Period	Not Required
Humidity	10 - 90% RH non-condensing
Operational vibration	10G, 10-40Hz Max, any direction
EMC	Meets all EC harmonized standards

### INPUT / OUTPUT SIGNALS

Detector In	BNC
Number of inputs	3
Trigger In	Rising or Falling Edge pulse Make/break
Camera Trigger	TTL positive pulse
Communication Interface	Data and command transfer via 1Gbps ethernet cable
Software	Custom software compatible with Microsoft Windows Operating Systems for control and data archiving in various file formats

### MECHANICAL

Dimensions mm (w/d/h)	650 x 230 x 310 (without tripod)
Weight	16kg / 35lbs (without camera and lens)
Mount	3/8-16 UNC Female

### CONTROL UNIT

System Clock	10MHz quartz crystal controlled
Trigger Jitter	<1us

### MIRROR

Type	Optical flat elliptical surface silvered
Size (HxW) mm	135 x 85

### INTEGRATED CAMERA OPTIONS

#### Tracker<sup>Lite</sup> -A

Photron FASTCAM  
Mini AX100 – 1024x1024 @ 4,000pps

#### Tracker<sup>Lite</sup> -V

Phantom VEO4 10L  
1280x800 @ 5,200pps

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com)

[info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## Lightweight Optical Motion Analysis System

**Compact Optical Tracking  
Pedestal mount**

**Multiple Tracking Modes**

**Up to 150KG Payload**

**Multiple Tracking  
and Payload Camera  
Options**



The Specialised Imaging LOMA system is the next generation of long range projectile tracking systems.

The system consists of an Azimuth module, Elevation module and payload platform, each constructed from stainless steel or aluminium to ensure rigidity and precision tracking.

The advanced servo control provided accurate ( $\pm 100$  micro-radians) movement of payloads up to 150Kg.

The LOMA system can be used in Manual, Slave or Optical tracking modes. IRIG-B timecode provided by existing range receiver or own receiver.

### FEATURES

- ☐  $\pm 100$  micro-radians angular accuracy
- ☐ Unlimited Azimuth rotation
- ☐ Elevation  $-20^\circ$  to  $+190^\circ$
- ☐ Turn and Plunge capable
- ☐ Ethernet communications
- ☐ Built in camera power, communications and trigger
- ☐ IR and Visible spectrum tracking options
- ☐ Customised Platforms

### OPERATING MODES

Manual	Joystick control for emplacement and maintenance purposes
Slave	External input of Longitude, Latitude & Elevation
Optical	Onboard camera tracking (Visible or IR)

### INPUT / OUTPUT SIGNALS

Command interface	RS232 (Controller to Pedestal)
Control interface	Gigabit Ethernet
Payload interface	Gigabit Ethernet, Coaxial.
Tracking interface	Video (tracking camera to Controller)

Azimuth/Elevation transfer RS232

Timecode interface	IRIG-B from existing generator
Software	Custom software compatible with Microsoft Windows Operating Systems for Control of Pedestal and Tracking cameras
Electrical Power	AC: 110-240V AC 50/60Hz or DC: 24-48V 30A

### OPERATING PARAMETERS

System Accuracy	$\pm 100$ micro-radians, $2\sigma$
Travel	Azimuth: Continuous (360°) Elevation: -20° to +190°
Torque (azimuth & elevation)	27Nm continuous (40Nm peak)
Nominal Payload	150Kg (at rated payload Mol)
Angular Velocity (nominal payload)	Azimuth: 90° / second Elevation: 90° / second
Angular Acceleration (nominal payload)	Azimuth: 90° / second <sup>2</sup> Elevation: 90° / second <sup>2</sup>
Non-Orthogonality	$\pm 10$ arc seconds
Drive System	Direct drive
Encoders	25-bit absolute, both axis

### ENVIRONMENTAL

Storage Temperature	-120°C to +70°C
Operating Temperature	-10C to +60°C
Humidity	20-80% RH non condensing
Vibration Shock	10-40 Hz Max. 10g in any direction
EMC	Meets all EC harmonised standards

### MECHANICAL

Dimension (w/d/h)	<b>Basic Pedestal:</b> 41cm x 39cm x 78.6cm (16.1" x 15.4" x 30.9") <b>Control Unit:</b> 72cm x 60cm x 50cm (28.3" x 23.6" x 19.7")
Weight	<b>Basic Pedestal:</b> 70Kg (154lbs) without base or Payload <b>Control Unit:</b> 35Kg (77lbs)
Transit Locks	Locking stow pins in Azimuth & Elevation

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## Multi-head Intensified Camera System

**Up to 8 camera heads per  
control module**

**1360 x 1024 pixel, 12-bit  
sensor resolution**

**Single or double  
image capture**

**5ns minimum  
shutter**



The Specialised Imaging CERBERUS camera system offers framing camera image capture performance with the addition of multiple camera control.

Each camera head can capture either one or two 1.4 MegaPixel 12-bit images with exposure times down to 5ns.

A maximum of 8 Control modules can be connected to allow a total of 64 cameras controlled from a single PC.

The CERBERUS system is flexible enough to allow multiple 3D/Stereoscopic image pairs or sequential images with a 5ns interframe time, equating to 200 Million Frames/second.

### FEATURES

- ☐ Control up to 64 camera heads
- ☐ Adjustable exposure down to 5ns
- ☐ Head to Head adjustable interframe time down to 5ns
- ☐ Nikon lens mount fitting
- ☐ Ethernet communications
- ☐ Compact and rugged design



### OPTICAL

Lenses	Nikon F-mount (ruggedized mounting system)
Shutter	Electro-mechanical
Distortion	Nominally zero
Intensity variation	Better than 5% across the image

### INTENSIFIER / SENSOR

#### CE-01 HEAD CE-02 HEAD

Image Sensor	ICX285AL (Intensified)	
Active CCD Pixel	1360 (H) x 1024 (V) pixels	
Pixel Size	6.45 µm (H) x 6.45 µm (V)	
Dynamic Range	12 bits	
Intensifier	8mm High resolution MCP Input window Fused Silica Output window Fibre Optic Photocathode S25, others on request Phosphor screen P43	
Gain	Up to 7,000	Up to 4,000
Dynamic resolution	>50lp/mm	>36lp/mm
Images	Single	Two (550ns interframe time)

### MECHANICAL

Dimension mm (w/d/h)	<b>Head (without lens)</b> 9.4cm x 21cm x 9.4cm (3.7" 8.2" x 3.7") <b>Controller</b> 19" rack mount 3U case
Weights	<b>Head</b> 3kg (6.6lbs) <b>Controller</b> 7kg (15.4lbs)
Head Mounting	3/8- 16 UNC Female in head base

### TIMING PARAMETERS

System Clock	200MHz quartz crystal control
Inherent Delay	500ns
Exposure Mode (each head)	Single exposure or multiple exposures (Max. 8) per head
Exposure Time	5ns – 10ms in 5ns steps
Interframe Time (head-to-head)	5ns – 20ms in 5ns steps
Delay to 1st exposure	500ns – 10ms in 5ns steps
Flash Outputs	5ns to 1ms in 5ns steps
Separation Time	30ns – 20ms in 5ns steps (multiple exposures on same channel)

### INPUT / OUTPUT SIGNALS

Trigger 1	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Trigger 2	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Timing Monitor Pulse	Pulse width (min. 5ns) and position user programmable TTL into 50Ω
Flash Trigger Outputs	Pulse width (min. 5ns) and position user programmable TTL into 50Ω
Remote Camera Interface	Data and command transfer via custom 10m cable.
Camera head control	Data and command transfer via 100Mbps Ethernet cable length 10m (standard), other lengths up to 100m available 100FX Fibre optic Ethernet link (up to 2Km) - optional
Software	Custom software compatible with Microsoft Windows Operating Systems for cameracontrol, image data archiving in various file formats.
Power Requirements	100-240V AC 2A, 50-60Hz

### ENVIRONMENTAL

Housing	19" Rack Mount 2U case
Storage temperature	-10°C to +50°C
Operating Temperature	-5°C to +40°C
Humidity	10 - 90% RH non condensing
Vibration shock	10 - 40 Hz Max. 10g in any direction
EMC	Meets all EC harmonized standards

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com)

[info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## High Resolution non-intensified camera

**Single image capture**

**Up to 11 MegaPixel  
12-bit images**

**Compact and rugged  
construction**



The Specialised Imaging T-Cam Camera offers single high resolution images for ambient or low light capture.

Comprehensive triggering adjustment and a wide range of output signals are controlled using the custom software package which includes measurement and image enhancement functions.

### FEATURES

- ☐ Fully adjustable exposure from 1 $\mu$ s to 10ms
- ☐ Adjustable output triggers
- ☐ Nikon lens mount fitting
- ☐ Gigabit ethernet communications

### OPTICAL

Lenses	Nikon F-mount
System Aperture	Limited by lens
Distortion	Nominally zero

### INTENSIFIER / SENSOR

	<b>T-Cam 43100</b>	<b>T-Cam 22100</b>
Image Sensor	KAI11002M (Non intensified)	KAI4021M (Non intensified)
Active CCD Pixel	4008 (H) x 2688 (V)	2048 (H) x 2048 (V)
Pixel Size	9µm (H) x 9µm (V)	7.4 µm (H) x 7.4 µm (V)
Dynamic Range	12 bits	12 bits

### MECHANICAL

Dimension mm (w/d/h)	17.8cm x 18.5cm x 19.7cm (7.0" x 7.3" x 7.8") without lenses
Mount	1/4-20 UNC Female (standard tripod)
Weight	6Kg (13.2lbs) without lens.

### INPUT / OUTPUT SIGNALS

Trigger 1	Electrical signal (BNC connector) Threshold variable from $\pm 25V$ Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Trigger 2	Electrical signal (BNC connector) Threshold variable from $\pm 25V$ Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Flash Trigger Output	Pulse width (min. 10ns) and position user programmable. TTL into 50Ω
Camera Control	Data and command transfer via Gigabit Ethernet. Cable length 10m (standard), other lengths up to 100m available 1000FX fibre optic Ethernet link (up to 2Km) - optional
Software	Custom software compatible with Microsoft Windows Operating Systems for camera control, image data archiving in various file formats.
Power Requirements	100-240V AC 2A, 50-60Hz

### ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10—90% RH non condensing
Vibration shock	10—40 Hz Max. 10g in any direction
EMC	Meets all EC harmonised standards

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## Compact image intensifier/converter

**Range of Intensifier tube options**

**Up to 40lp/mm resolution**

**Gate times down to 50ns**

**Comprehensive triggering  
for High-speed video  
and Standard  
cameras**



The Specialised Imaging SIL3 image intensifier is the latest generation of intensifier that can be synchronised to almost any high-speed video, conventional video or stills camera.

Maximum equivalent frame rates ranging from 100,000fps to 1,000,000fps and minimum gate (exposure) time of 50ns.

The range of Intensifier tube options include output window diameters, Phosphor decay times and Wavelength responses (inc. UV) provide a unique enhancement to any camera.

Intuitive custom control software allows gain, gate/exposure and delay adjustment via Ethernet.

### FEATURES

- ☐ 25mm and 40mm Intensifier tube options
- ☐ Fully adjustable gating time down to 50ns
- ☐ Gain adjustment up to 500,000X
- ☐ Nikon lens mount fitting
- ☐ Optical camera coupling
- ☐ Ethernet communications
- ☐ UV wavelength conversion option

## INTENSIFIER /SENSOR

	SIL40HG50	SIL40NG50	SIL25HG50	SIL25NG50	SIL25HG50-X	SIL25HG50-D
Front Window	Fused Silica	Fused Silica	Fused Silica	Fused Silica	Fused Silica	Quartz
Output Window	Glass	Glass	Glass	Glass	Glass	Glass
Photocathode	S20	S20	S20	S20	S20	UV Enhanced S20
Output Phosphor	FS/10µs decay	FS/10µs decay	FS/10µs decay	FS/10µs decay	P46/300ns decay	P46/300ns decay
Input Area (diameter)	40mm	40mm	25mm	25mm	25mm	25mm
Spectral response (minimum)	15mA/W @ 214nm	17mA/W @ 214nm	15mA/W @ 214nm	20mA/W @ 214nm	20mA/W @ 214nm	20mA/W @ 220nm
	55mA/W @ 450 nm	55mA/W @ 450 nm	48mA/W @ 450 nm	53mA/W @ 450nm	53mA/W @ 450nm	53mA/W @ 450nm
	17mA/W @ 800 nm	10mA/W @ 800 nm	17mA/W @ 800 nm	17mA/W @ 800 nm	17mA/W @ 800 nm	17mA/W @ 800 nm
Gain (W/W @ 500nm)	100,000	10,000	100,000	10,000	100,000	500,000
Limiting resolution (Typical)	22 lp/mm	30 lp/mm	27 lp/mm	40 lp/mm	27 lp/mm	28 lp/mm
Frame Rate (Max. Fps)	100,000	100,000	100,000	100,000	1,000,000	1,000,000

## OPTICAL

Front Lens mount	F-mount
Optical Coupling	F-mount

## MECHANICAL

Dimension mm (w/d/h)	17.8cm x 18.5cm x 19.7cm (7.0" x 7.3" x 7.8") without lenses
Mount	3/8- 16UNC Female (standard tripod)
Weight	4.8Kg (10.6lbs) without lenses
Lens/Camera Support	Option available

## ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5° to +40°C
Humidity	10 - 90% RH non condensing
Vibration shock	10 - 40Hz Max. 10g in any direction
EMC	Meets all EC harmonized standards

## TIMING PARAMETERS

System Clock	200MHz, quartz crystal controlled.
Inherent Delay (5ns jitter)	60ns
Exposure Time	50ns—DC in 5ns steps independently variable
Inter Exposure Time	50ns—25ms in 5ns steps independently variable
Number of exposures	Up to 64 in pulsed mode. Unlimited in REP (synchronous) mode

## INPUT / OUTPUT SIGNALS

Triggers	Electrical signal (BNC connector) Maximum Input level 50V Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Aux Outputs	Pulse width and position user programmable (min. 5ns) TTL into 50Ω
Control Interface	Remote control via Standard 100Mbps Ethernet
Inhibit Input	5V TTL (user input brightness protection)
Software	Custom software compatible with Microsoft Windows Operating Systems for control
Power Requirements	100-240V AC 2A, 50-60Hz


**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.





## High Resolution Dual Image Intensified camera

**Single or Double image capture**

**Up to 10.7 MegaPixel  
12-bit images**

**Lightweight and  
Rugged construction**



The Specialised Imaging SIR3 Framing Camera offers up to 2 high resolution images, 100µs apart. Fully flexible intensified CCD sensor provides control over interframe time, gain and exposure.

Comprehensive triggering adjustment and a wide range of output signals are controlled using the custom software package which also includes measurement and image enhancement functions.

### FEATURES

- ☐ Fully adjustable interframe time to 100µs
- ☐ Fully adjustable exposure down to 10ns
- ☐ Gain adjustment up to 10,000X
- ☐ Adjustable output triggers
- ☐ Nikon lens mount fitting
- ☐ Gigabit ethernet communications

### OPTICAL

Lenses	Nikon F-mount (ruggedized mounting system)
System Aperture	f 2
Shutter	Electro-mechanical
Distortion	Nominally zero
Coupling	CCD to MCP via FO
Vignetting	<3%
Intensity variation	Better than 5% across the image
Optical Viewfinder	Optional

### INTENSIFIER / SENSOR

	<b>SIR3-18D</b>	<b>SIR3-25D</b>	<b>SIR3-40D</b>
Image Sensor	ICX285AL	KAI4021M	KAI11002M
Active CCD Pixel	1360 (H) x 1024 (V)	2048 (H) x 2048 (V)	4008 (H) x 2688 (V)
Pixel Size	6.45 µm (H) x 6.45 µm (V)	7.4 µm (H) x 7.4 µm (V)	9 µm (H) x 9 µm (V)
Dynamic Range	12 bits	12 bits	12 bits
Intensifier diameter	18mm MCP	25mm MCP	40mm MCP
Photocathode	All models: S25		
Phosphor / decay	P46/300ns	FS/10µs	FS/10µs
Input / Output windows	All models: Glass / Fibre		
Gain	Variable up to 10,000 all models		

### MECHANICAL

Dimension mm (w/d/h)	17.0cm x 48.5cm x 19.3cm (without lens)
Mount	1/4 - 20 UNC and 3/8 - 16 UNC female
Weight	15Kg (33lbs) without lens

### ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10—90% RH non condensing
Vibration shock	10—40 Hz Max. 10g in any direction
EMC	Meets all EC harmonised standards

### TIMING PARAMETERS

System Clock	200MHz quartz crystal controlled
Inherent Delay	<130ns
Imaging Mode	Single or Double image
Exposure Modes (each image)	Single exposure or multiple exposures (Max. 16 - subject to imaging conditions).
Exposure Times	10ns – 10ms in 5ns steps independently variable
Delay to 2nd exposure	100µs – 10mS in 5ns steps.
Flash output	20ns to 1ms in 5ns steps independently variable
Separation	30ns to 20ms in 5ns steps independently variable

### INPUT / OUTPUT SIGNALS

Trigger 1	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Trigger 2	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination
Flash Trigger Output	Pulse width (min. 10ns) and position user programmable. TTL into 50Ω
Camera Control	Data and command transfer via Gigabit Ethernet Cable length 100m (standard) 1000FX fibre optic Ethernet link (up to 2Km) - optional
Software	Custom software compatible with Microsoft Windows Operating Systems for camera control, image data archiving in various file formats.
Power Requirements	100-240V AC 2A, 50-60Hz

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**


**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com)

[info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



## Ballistic Velocity Measurement System

**3,330 measurements /  
second**

**Up to 15,000  
measurements storage**

**5ns Resolution**

**3 independent input  
channels**

**4 independent output  
channels**



The Specialised Imaging VT system performs time of flight and/or velocity measurements using 3 input channels.

Capture rates up to 3,330 measurements per second (200,000 rounds per minute) can be measured and up to 15,000 measurements can be stored locally in the Unit head.

The four independent Output channels can be used as a trigger adjustment unit for ballistic events to ensure correct triggering of framing, ultra high-speed and high speed video cameras.

### FEATURES

- ☐ Velocity measurement
- ☐ Time of flight measurement
- ☐ Ethernet communication

## OPERATING PARAMETERS

Measurements	Time of flight or Velocity
Max. rate of capture/fire	200,000 measurements/minute (3,330 Hz)
Max. storage	15,000 measurements (velocity or time of flight)

## INPUT / OUTPUT SIGNALS

Input/trigger channels	3 (independent) Electrical signal (BNC connector) Threshold variable from $\pm 25V$ Positive or Negative polarity, Make/Break 50 $\Omega$ or 1K $\Omega$ termination
Output channels	4 (independent) +ve TTL (BNC connector) Pulse width user programmable 50 $\Omega$ termination
Control interface	100Mbps Ethernet
Software	Custom software compatible with Microsoft Windows Operating Systems for Control and data archiving.
Electrical power	AC: 110-240V AC 50/60Hz or DC: 24V DC – Battery power

## TIMING PARAMETERS

System Clock	200MHz, quartz crystal controlled
Inherent Delay	50ns
Range	5ns to 1s in 5ns steps independently variable
Jitter (trigger)	$\leq 5ns$
Jitter (channel to channel)	0ns

## MECHANICAL

Dimension mm (w/d/h) Inc. handles	23cm x 23 cm x 17cm (9.1" x 9.1" x 6.7")
IP rating	IP65
Weight	4Kg (8.8lbs)

## ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating Temperature	-5°C to +40°C
Humidity	10 - 90% RH non condensing
Vibration shock	10 - 40 Hz Max. 10g in any direction
EMC	Meets all EC harmonized standards


**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com)

[info@specialised-imaging.com](mailto:info@specialised-imaging.com)



## Laser illumination system

Up to 400W lighting power

640nm wavelength

Flexible triggering/  
synchronisation



The Specialised Imaging LUX640 laser illumination system provides up to 400W of lighting power at pulse frequencies up to 10MHz or single pulses up to 30 $\mu$ S.

Simple triggering allows the SI-LUX640 to interface with most high-speed cameras ranging from High-speed video to Ultra highspeed framing cameras.

The 2m laser output light guide includes user interchangeable low coherence beam expanders.

### FEATURES

- ☐ Low coherence
- ☐ Pulse width from 10ns – 30 $\mu$ s
- ☐ Pulse frequency range from single to 10MHz
- ☐ Compact design



### OPTICAL

Beam Expanders	φ25mm and φ50mm versions
Wavelength	640nm ± 6nm
Output power	200W (-10/ + 30%) or 400W (-10/ + 20%)
Power drop	~ 0.2% / μS for pulses less than 5μS

### INPUT / OUTPUT SIGNALS

Sync. Input	5V TTL (BNC socket connector) (laser pulse duration = duration of +5V state)
Indicators	Green LED – Laser is powered & ready
Software	Custom software compatible with Microsoft Windows Operating Systems for Control.
Electrical Power	100-240V AC 2A 50-60Hz

### HEALTH & SAFETY

Laser rating	Class 3b
Safety features	Key operated master control power on/off Connector for remote Interlock
Visual indicator	Green LED = Laser is powered & ready Red LED = interlock indicator Health & Safety on laser body

### ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10—90% RH non condensing
Vibration shock	10—40 Hz Max. 10g in any direction
EMC	Meets all EC harmonised standards

### TIMING PARAMETERS

Max pulse rate	10MHz
Min. single pulse duration	~10ns
Max. single pulse duration	30μS (max. power drop 20%)
Rise time	~10ns (10%...90%)
Fall time	~5ns
Delay	70ns between input to start of light pulse (incl. control cable delay)
Jitter	<5ns
Operating limits	Max. 0.03% duty cycle for unlimited operation Max. 100% duty cycle for 30μS laser operation

### MECHANICAL


Dimension (w/d/h)	<b>Laser:</b> 6.2cm x 15cm x 3.6cm (6" x 2.4" x 1.4") <b>Master control box:</b> 12cm x 6cm x 12cm (4.7" x 2.4" x 4.7")
Weight	<b>Laser:</b> 0.5Kg (1.1lbs) <b>Master control box:</b> 0.5Kg (1.1lbs)
Light guide	2m
Control cable	2m between control/safety box & laser

**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)



# SI-OT3



## Optical Trigger Unit

**Rugged Design**

**Battery Powered**

**Shadow and IR Flash Detector**

**High / Low pass filtering**



The Specialised Imaging OT3 provides a reliable optical trigger for either projectile "shadow" detection or IR flash detection

Battery powered, with a rugged enclosure allow the OT3 to be used outside in all weathers and independent of mains power.

Battery unit includes mains charger and can power the OT3 for up to 8 hours.

### FEATURES

- ☐ Small and lightweight
- ☐ Battery powered for up to 8 hours
- ☐ Nikon lens mount fitting
- ☐ User adjustable sensitivity

### OPTICAL

Lenses	Nikon F-Mount
Alignment	Optical viewport

### INTENSIFIER / SENSOR

Sensor	Multi-Segment Photodiode array. 300nm – 700nm range (Non-intensified)
--------	--

### INPUT / OUTPUT SIGNALS

Output	Positive 5V TTL (BNC socket connector) 50Ω termination
Trigger indicator	LED
Software	Custom software compatible with Microsoft Windows Operating Systems for Control and data archiving.
Electrical Power	DC: 18-34V DC Battery powered: Sealed lead acid Built in battery charger

### MECHANICAL

Dimension mm (w/d/h)	<b>Sensor head (without lens)</b> 9cm x 17cm x 9cm (3.5" x 6.7" x 3.5") <b>Battery</b> 15cm x 9.5cm x 17cm (6" x 3.7" x 6.7")
Weights	<b>Sensor head</b> 1Kg (2.2lbs) <b>Battery</b> 4Kg (8.8lbs)
Sensor head mounting	1/4 -20 UNC UNC Tripod Female in base & top

### ENVIRONMENTAL

Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10—90% RH non condensing
Vibration shock	10—40 Hz Max. 10g in any direction
EMC	Meets all EC harmonised standards


**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com) [info@specialised-imaging.com](mailto:info@specialised-imaging.com)


As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.



# SI-AT1



## Acoustic trigger Unit



**Built in battery power**

**Compact size**

**User selectable sensitivity settings**

**TTL trigger output**

The Specialised Imaging SI-AT1 provides a reliable acoustic trigger for blast or shockwave detection.

Built in battery power within a compact, rugged enclosure allows the AT1 to be used outside in all weathers and independent of mains power for up to 8 hours.

### FEATURES

- ☐ Small and rugged
- ☐ Battery powered for up to 8 hours
- ☐ User adjustable sensitivity
- ☐ Standard tripod mounting threads

### SENSOR

Sensor	Piezoelectric acoustic sensor
--------	-------------------------------

### INPUT / OUTPUT SIGNALS

Output	Positive 5V TTL (BNC socket connector) 50Ω termination
Pulse width	10mS
Trigger indicator	LED
Electrical power	Battery powered Built in battery charger 12VDC - 500mA (from charger)
Lower power indicator	LED

### MECHANICAL

Dimension mm (w/d/h)	<b>Unit</b> 152mm x 88mm (5.9" x 3.46") <b>Tripod mount block</b> 16mm high (0.62")
Weight	1Kg (2.2lbs)
Tripod mount	1 x 1/4-20UNC 1 x 3/8-16UNC

### ENVIRONMENTAL

Storage temperature	-20°C to +50°C
Operating temperature	-15°C to +50°C
Charging temperature	0°C to +40°C
Humidity	10—90% RH non condensing
Vibration shock	10—40 Hz Max. 10g in any direction
EMC	Meets all EC harmonised standards


**UK** (Head Office / Factory)  
6 Harvington Park,  
Pitstone Green Business Park  
Pitstone. LU7 9GX England  
**Tel +44 (0) 1442 827728**

**USA**  
Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA  
**Tel +1 951-296-6406**

**GERMANY**  
Hauptstr. 10,  
82275 Emmering  
Germany  
**Tel +49 8141 666 89 50**

[specialised-imaging.com](http://specialised-imaging.com)    [info@specialised-imaging.com](mailto:info@specialised-imaging.com)

As part of our on-going commitment to improvement we reserve the right to alter specifications, designs or figures, without prior notice. All dimensions and weights are approximate.





# AD-500



## Single or multiple head high intensity flash system



**500J Flash head with 2ms duration**

**Up to four independent flash heads**

**40s recycle time**

**Standard trigger**

The Specialised Imaging AD500 Flash system offers the flexibility of four controllable high intensity flash light for used in scientific and industrial environments.

### CONTROL UNIT

Signal Input	Four independent Channels
Trigger Mode	1. Independently 2. Simultaneous - all four channels triggered via channel 1
Trigger Source	Short Circuit 5-100V positive edge
Input Impedance	50Ω per channel
Mains Input	IEC socket
Input Voltage	90-240V 50-50Hz AC
Dimensions mm	(LxWxH) 220mm x 110mm x 128mm
Weight	3.7kg
Lamp cable length	2.5m


### FLASH HEAD

Light Duration (Typ)	2ms measured to 50% of peak output
Stored Charge (max)	500J
Charge Voltage	340V
Light Source	U-Shape Xenon flashtube
Rise Time	50μs
Delay (typ)	30μs
Recycling Time (typ)	40 seconds
Dimensions	(LxDia) 270mm x 170mm
Weight (Kg)	5.25kg

# MSFH-370



**High intensity flash with dedicated fibre optic output**



**370J Flash head with 750µs duration**

**Up to four independent flash heads**

**40s recycle time**

**Standard trigger**

The Specialised Imaging MSFH-370 Flash system offers the flexibility of four controllable high intensity flash light for used in scientific and industrial environments. 5mm diameter FO output.

## CONTROL UNIT

Signal Input	Four independent Channels
Trigger Mode	1. Independently 2. Simultaneous - all four channels triggered via channel 1
Trigger Source	Short Circuit 5-100V positive edge
Input Impedance	50Ω per channel
Mains Input	IEC socket
Input Voltage	90-240V 50-50Hz AC
Dimensions mm	(LxWxH) 220 x 120 x 90 mm
Weight	4.5 kg

## FLASH HEAD

Light Duration (Typ)	750µs measured to 50% of peak output
Stored Charge (max)	370J
Charge Voltage	340V
Light Source	Linear spark source flashtube
Rise Time	50µs
Delay (typ)	30µs
Recycling Time (typ)	40 seconds
Dimensions	(LxDia) 270 x 170 mm
Weight (Kg)	5.25 kg (without handle, legs, FO mounting)



Ballistics  
Detonics  
Plasma  
Impact studies  
Combustion research  
Spray and particle analysis  
Medical testing and research  
Low light machine vision system  
Nanotechnology and micro-machines  
Elasticity, crack propagation and shock resistance

**[www.specialised-imaging.com](http://www.specialised-imaging.com)**

**[info@specialised-imaging.com](mailto:info@specialised-imaging.com)**



FM 87429 ISO9001:2015



THE QUEEN'S AWARDS  
FOR ENTERPRISE:  
INNOVATION  
2016

**UK (HEAD OFFICE / FACTORY)**

6 Harvington Park,  
Pitstone Green Business Park,  
Pitstone, LU7 9GX. United Kingdom

**Tel +44 (0) 1442 827728**  
**Fax +44 (0) 1296 668098**

**USA**

Specialised Imaging Inc.  
40935 County Center Dr. Suite D  
Temecula, CA 92591, USA

**Tel +1 951-296-6406**

**GERMANY**

Hauptstr. 10,  
82275 Emmering  
Germany

**Tel +49 8141 666 89 50**  
**Fax +49 8141 666 89 33**