

Digital Microscope



# Powerful Analysis, Dynamic Imaging







# Intelligent Innovation

# Macro to Micro

# Large selection of lenses to find the best magnification, resolution, and working distance for your sample


- 23X to 8220X magnification
- Accommodates larger samples with an up to 66 mm (2.6 in.) working distance

### Free angle observation system

- Oblique observation (± 90°)
- XY motorized stage with rotation (± 90°)



17 available lenses



# Be Confident in Your Results with Guaranteed Accuracy and Precision

#### Precise measurements

• Eliminate the image conversion effect for accurate measurements with a telecentric optical system

### Accuracy and repeatability guaranteed at all


### magnifications

· Acquire reliable measurement results



Reliable, accurate measurements

\*To issue certificates, calibration work must be undertaken by Olympus's dedicated service technician.



some manual

XY mo stage rotation

Obli

(± §

obser

ammun



Remote console makes the system easy to control


## **Multiple Observations with a Single Click**

#### Instant switching saves time

- Quickly change objectives
- Switch observation methods by pushing a button

#### All observation methods available at all magnifications

• Six observation functions are built into the optical head



Switching between observation methods is as push of a button


# Advanced Measurements Are Fast and Easy to Obtain

# Improved analysis functions make the DSX1000 powerful and versatile

Measure profiles, inter-surface steps, surface roughness, and more

#### Faster analyses with advanced easy-to-use functions

· Automated functions ease your workflow



Powerful analysis software

# Specifications Main Unit Specifications

		DSX10-SZH	DSX10-UZH	
Optical system	Optical system	Telecentric optical system		
	Zoom ratio	10X (mortorized)		
	Zoom magnification method	Motorized		
	Calibration	Automatic		
	Lens attachment	Quick-switch, coded lens attachments automatically update magnification and visual field information		
	Maximum total magnification (on a 27-inch monitor monitor)	8220X		
	Working distance (W.D)	66.1 mm – 0.35 mm		
	Accuracy and repeatability (X-Y plane)	Accuracy of magnification: 3% <sup>*1</sup> Repeatability of magnification: ± 2% = 3 <sub>0</sub> n - 1		
	Repeatability (Z axis)	Repeatability (height): $\sigma n - 1 \le 1 \mu m^{*2}$		
	Image sensor	1 / 1.2 inch, 2.35 million pixel color CMOS		
	Cooling	Peltier cooling		
Camara	Frame rate	60 fps (maximum)		
Camera	Normal	1,200 × 1,200 (1:1) / 1,600 × 1,200 (4:3)		
	Fine	Not available	1,200 × 1,200 (1:1) / 1,600 × 1,200 (4:3)	
	Super fine	Not available	3,600 × 3,600 (1:1) / 4,800 × 3,600 (4:3)	
	Color light souce	LED		
Illumination	Lifetime	60,000 h (design value)		
	BF (brightfield)	Standard		
	OBQ (oblique)	Standard		
	DF (darkfield)	Standard LED ring divided into four divisions		
Observation	MIX (brightfield+darkfield)	Standard Simultaneous observation of BF + DF		
	PO (polarization)	Standard		
	DIC (differential interference)	Not available	Standard	
	Contrast up	Standard		
	Depth of focus up function	Not available	Standard	
	Transmitted lighting	Standard*3		
Focus	Focusing	Motorized		
	Stroke	101 mm (motorized)		

\*1 Calibration by Olympus or dealer service technician necessary. \*2 When used 20X or higher objective. \*3 The optional DSX10-ILT is required.

Objective		DXS10-SXLOB	DSX10-XLOB	UIS 2
Objective lens	Maximum sample height	50 mm (2 in.)	115 mm (4.5 in.)	145 mm (5.7 in.)
	Maximum sample height (free angle observation)	50 mm (2 in.)		
	Parfocal distance	140 mm (5.5 in.)	75 mm (3 in.)	45 mm (1.8 in.)
	Lens attachment	Integrated with lens	Available	
	Total magnification (on a 27-inch monitor)	23–1644x	49–6570x	26* <sup>4</sup> -8220x
	Actual F.O.V. (µm)	19,200 μm – 270 μm	9,100 µm – 70 µm	17,100 μm – 50 μm
Adaptor	Diffusion adaptor (optional)	Available	Not available	
	Eliminate reflection adaptor (optional)	Available	Not available	
Lens attachment	Number of objectives that can be attached	Up to 1 piece (attachment is integrated with lens)	Up to 2 pieces	
Objective lens case		Three lens attachments can be stored		

\*4 Total magnification when using MPLFLN1.25X

Stage		DSX10-RMTS	DSX10-MTS	DSX10-MTS	
Stage	XY stage: motorized / manual	Motorized (with rotation function)	Motorized	Motorized	
	XY stroke	Stroke priority mode : 100 mm × 100 r Rotation priority mode : 50 mm × 50 n		ım	100 mm × 105 mm
	Rotation angle	Stroke priority mode : ±20° Rotation priority mode : ±90°		Not available	
	Display rotation angle	GUI		Not available	
	Load-resistance	5	i kg (11 lb)		1 kg (2.2 lb)
Frame	Upright frame	Tilt frame		Display	у
Z-axis stroke	troke 50 mm (manual)		Size		27-inch flat panel display
Tilt observetio	n Not available	±90°	Resolution		1,920 (H) × 1,080 (V)
Tilt angle displa	ay Not available	GUI			
Tilt angle meth	od Not available	Manual, fix / release handle			

System Total	Upright frame system	Tilt frame system	
Weight (frame, head, motorized stage, display, and console)	43.7 kg (96.3 lb)	46.7 kg (103 lb)	
Power consumption	100 V-120 V/220 V-240 V-1.1/0.54 A-50/60 Hz		

#### • EVIDENT CORPORATION is ISO14001 certified.

• EVIDENT CORPORATION is ISO9001 certified.

All company and product names are registered trademarks and/or trademarks of their respective owners. Olympus and the Olympus logo are trademarks of Olympus Corporation or its subsidiaries.
Performance characteristics and other values described in this brochure are based on Olympus evaluations as of October 2021, and are subject to change without notice.
The information including guaranteed accuracy in this brochure is based on the condition set by Olympus. For details, refer to the Instruction Manual. Images on the PC monitors are simulated.
Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

## EvidentScientific.com



EVIDENT CORPORATION Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0910, Japan

